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
The screenshot shows the 2100 Expert software interface with several key sections highlighted by red boxes and numbered instructions:

- Instrument Identification:** A photo of the DNA Chip is shown on the left. To its right, fields for Name, Serial#, Cartridge, Vendor, and Product ID are visible. A red box highlights the 'Simulation Mode' firmware status with the instruction: "3. 确认芯片已被识别。"
- COM Port:** A dropdown menu is set to '1', with a red box and instruction: "1. 确保选定正确端口。"
- Assay Selection:** The 'Assays' button is highlighted with a red box and instruction: "2. 点击选择试验类型。"
- Start Button:** The 'Start' button is highlighted with a red box and instruction: "7. 点击开始运行。"
- Assay File:** The file path 'C:\...sis\Demo DNA 7500 Series II.xs...' is highlighted with a red box.
- Destination and Data Acquisition Parameters:** The 'Default' radio button is selected. The file path 'C:\...ilent\2100 bioanalyzer\2100 expert\Data\2009-05-20' is highlighted with a red box and instruction: "5. 设置数据保存路径以及文件名前缀。"
- File Prefix:** The text '2100 expert' is entered in the field, highlighted with a red box.
- Run sample:** The range '1 to 12' is set in the spinner controls, highlighted with a red box.
- Assay Details:** The 'Assay Class' is set to 'DNA 7500', highlighted with a red box.
- Start Run Checklist:** A list of seven items is shown, all with green checkmarks. A red box highlights the entire checklist with the instruction: "4. 检查所有前提条件已经满足。"
- Chip Summary Table:** A table with columns for Sample Name, Sample Comment, Rest. Digest, and Observation. The first row is 'ladder 1'. A red box highlights the table with the instruction: "6. 填写样品名，也可直接开始，数据获取后可以更改样品名。"

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The screenshot shows the 'Electropherogram - Ladder 1' software interface. The main window displays two electropherogram plots. A 'Result Flagging' dialog box is open, showing a table of fragments to search for and options to specify results. A 'Print' dialog box is also open, showing options to print various items and save to file. Chinese annotations are present throughout the image, providing instructions on how to use the software.

Annotations in the image include:

- 1. 点击打开数据 (Click to open data)
- 2. 设定结果标记规则 (Set result marking rules)
- 3. 点击生成报告 (Click to generate report)
- 4. 应用设定的规则 (Apply the set rules)
- 5. 点击Save生成PDF报告 (Click Save to generate PDF report)

The 'Result Flagging' dialog box contains the following information:

1. Specify the fragments to search for

Index	Fragment Size [bp]
1	100
2	200

2. Specify options

in a tolerance of 5 percent

3. Specify results

The samples that contain the specified fragments

Shall be colored: Change...

Shall be labelled: Sample contains 100bp and 200bp fragments

All other samples

Shall be colored: Change...

Shall be labelled: All Other Samples

The 'Print' dialog box contains the following information:

Print Item

Run Summary Electropherogram Calibration Curve Run Logbook

Assay Details Gel Like Standard Curve

Result Tables

Wells

All Wells

Wells [1,4-9,12]

Options

2 per page

Exclude Marker

Include Ladder

Save To File

PDF File Path: C:\...2100 expert\Data\Demo DNA 1000 Series II.pdf

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